

# Product Data Sheet

## Anti-CCL8 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1065-C	Mouse	Human	SE
<b>Description</b>	Mouse monoclonal capture antibody to CCL8		
<b>Immunogen</b>	Recombinant protein of human CCL8		
<b>Purification</b>	The antibody was purified by Protein A/G affinity chromatography.		
<b>Specificity</b>	Recognizes native and recombinant human CCL8 protein		
<b>Clonality</b>	Monoclonal		
<b>Conjugation</b>			
<b>Form</b>	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
<b>Dilution</b>			
<b>Gene Symbol</b>	CCL8		
<b>Alternative Names</b>	MCP2; SCYA10; SCYA8; C-C motif chemokine 8; HC14; Monocyte chemoattractant protein 2; Monocyte chemotactic protein 2; MCP-2; Small-inducible cytokine A8		
<b>Entrez Gene</b>	6355		
<b>SwissProt</b>	P80075		
<b>Directions for Use</b>	Recommended pair for sandwich immunoassay - Detection: CAP1065-D		
<b>Storage/Stability</b>	Shipped at 4° C. Upon delivery aliquot and store at -20° C for 1 year. Avoid freeze/thaw cycles.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com