

# Product Data Sheet

## CD158b1 Blocking Peptide

Catalog #	Source	Reactivity	Applications
CBP1653	Synthetic	H, M, R	BL
<b>Description</b>	The peptide is used to block Anti-CD158b1 Antibody (#CPA1653) reactivity.		
<b>Form</b>	Lyophilized powder		
<b>Gene Symbol</b>	KIR2DL2		
<b>Alternative Names</b>	CD158B1; NKAT6; Killer cell immunoglobulin-like receptor 2DL2; CD158 antigen-like family member B1; MHC class I NK cell receptor; Natural killer-associated transcript 6; NKAT-6; p58 natural killer cell receptor clone CL-43; p58 NK receptor CL-43; CD158b1		
<b>Entrez Gene</b>	3803 (Human)		
<b>SwissProt</b>	Q9UM54 (Human); Q64331 (Mouse)		
<b>Purity</b>	>85%		
<b>Quality Control</b>	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.		
<b>Directions for Use</b>	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com