

# Product Data Sheet

## CSGALNACT1 Blocking Peptide

| Catalog #                 | Source                                                                                                                                                                            | Reactivity | Applications |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CBP2715                   | Synthetic                                                                                                                                                                         | H          | BL           |
| <b>Description</b>        | The peptide is used to block Anti-CSGALNACT1 Antibody (#CPA2715) reactivity.                                                                                                      |            |              |
| <b>Form</b>               | Lyophilized powder                                                                                                                                                                |            |              |
| <b>Gene Symbol</b>        | CSGALNACT1                                                                                                                                                                        |            |              |
| <b>Alternative Names</b>  | CHGN; GALNACT1; Chondroitin sulfate N-acetylgalactosaminyltransferase 1; CsGalNAcT-1; Chondroitin beta-1,4-N-acetylgalactosaminyltransferase 1; Beta4GalNAcT-1                    |            |              |
| <b>Entrez Gene</b>        | 55790 (Human)                                                                                                                                                                     |            |              |
| <b>SwissProt</b>          | Q9P218 (Human)                                                                                                                                                                    |            |              |
| <b>Purity</b>             | >85%                                                                                                                                                                              |            |              |
| <b>Quality Control</b>    | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.                                                                                         |            |              |
| <b>Directions for Use</b> | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Store at -20°C for one year.                                                                                                                                      |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com