

# Product Data Sheet

## RagA/B Blocking Peptide

Catalog #	Source	Reactivity	Applications
CBP4511	Synthetic	H, M, R, B, Mk	BL
<b>Description</b>	The peptide is used to block Anti-RagA/B Antibody (#CPA4511) reactivity.		
<b>Form</b>	Lyophilized powder		
<b>Gene Symbol</b>	RRAGA; RRAGB		
<b>Alternative Names</b>	RRAGA; Ras-related GTP-binding protein A; Rag A; RagA; Adenovirus E3 14.7 kDa-interacting protein 1; FIP-1; RRAGB; Ras-related GTP-binding protein B; Rag B; RagB		
<b>Entrez Gene</b>	10670, 10325 (Human); 68441, 245670 (Mouse); 117044, 100909655 (Rat)		
<b>SwissProt</b>	Q7L523, Q5VZM2 (Human); Q80X95, Q6NTA4 (Mouse); Q63486, Q63487 (Rat)		
<b>Purity</b>	>85%		
<b>Quality Control</b>	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.		
<b>Directions for Use</b>	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com