

Fibrinopeptide B, Human

Catalog #	Source	Reactivity	Applications
CCP1205	Synthetic		

Description Peptide to Fibrinopeptide B, Human

Biological Description

Form Lyophilized powder

CAS Number 36204-23-6

Molecular Formula C₆₆H₉₃N₁₉O₂₅


Molecular Weight 1552.59

Purity > 95%

Chemical Structure pGlu - Gly - Val - Asn - Asp - Asn - Glu - Glu - Gly - Phe - Phe - Ser - Ala - Arg - OH

Storage/Stability Shipped at 4°C. Store at -20°C for one year.

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com