

GLP-2, Human

Catalog #	Source	Reactivity	Applications
CCP1240	Synthetic		
Description	Peptide to GLP-2, Human		
Biological Description			
Form	Lyophilized powder		
CAS Number	99120-49-7		
Molecular Formula	C171H266N48O56S		
Molecular Weight	3922.38		
Purity	> 95%		
Chemical Structure	H - His - Ala - Asp - Gly - Ser - Phe - Ser - Asp - Glu - Met - Asn - Thr - Ile - Leu - Asp - Asn - Leu - Ala - Ala - Arg - Asp - Phe - Ile - Asn - Trp - Leu - Ile - Gln - Thr - Lys - Ile - Thr - Asp - Arg - OH		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com