


[Val438]-Tyrosinase (432-444) (Human)

Catalog #	Source	Reactivity	Applications
CCP1685	Synthetic		
Description	Peptide to [Val438]-Tyrosinase (432-444) (Human)		
Biological Description			
Form	Lyophilized powder		
CAS Number	320341-56-8		
Molecular Formula	C ₆₅ H ₉₃ N ₁₅ O ₂₆		
Molecular Weight	1500.55		
Purity	> 95%		
Chemical Structure	Ser - Tyr - Leu - Gln - Asp - Ser - Val - Pro - Asp - Ser - Phe - Gln - Asp		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com