


[Des-Ser1]-Cerebellin

Catalog #	Source	Reactivity	Applications
CCP1687	Synthetic		
Description		Peptide to [Des-Ser1]-Cerebellin	
Biological Description			
Form		Lyophilized powder	
CAS Number		94245-80-4	
Molecular Formula		C66H108N22O21	
Molecular Weight		1545.73	
Purity		> 95%	
Chemical Structure		Gly - Ser - Ala - Lys - Val - Ala - Phe - Ser - Ala - Ile - Arg - Ser - Thr - Asn - His	
Storage/Stability		Shipped at 4°C. Store at -20°C for one year.	
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com