

Product Data Sheet

[D-Arg1, D-Phe5, D-Trp7, 11]-Substance P

Catalog #	Source	Reactivity	Applications
CCP1761	Synthetic		
Description		Peptide to [D-Arg1, D-Phe5, D-Trp7, 11]-Substance P	
Form		Lyophilized powder	
CAS Number			
Molecular Formula		C75H101N19O12	
Molecular Weight		1460.76	
Purity		> 95%	
Chemical Structure		D - Arg - Pro - Lys - Pro - D - Phe - Gln - D - Trp - Phe - Gly - Leu - D - Trp - NH ₂	
Storage/Stability		Shipped at 4°C. Store at -20°C for one year.	
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com