

Product Data Sheet

[Des-Arg30, Des-Pro31]-Dendroaspis Natriuretic Peptide

Catalog #	Source	Reactivity	Applications
CCP1763	Synthetic		
Description	Peptide to [Des-Arg30, Des-Pro31]-Dendroaspis Natriuretic Peptide		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C169H263N51O54S2		
Molecular Weight	3937.42		
Purity	> 95%		
Chemical Structure	Glu - Val - Lys - Tyr - Asp - Pro - Cys - Phe - Gly - His - Lys - Ile - Asp - Arg - Ile - Asn - His - Val - Ser - Asn - Leu - Gly - Cys - Pro - Ser - Leu - Arg - Asp - Pro - Asn - Ala - Pro - Ser - Thr - Ser - Ala (Disulfide bridge Cys7 - Cys23)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com