

Product Data Sheet

[Gly1, Ser3, 22, Gln4, 34, Thr6, Ala19, Tyr21, Ala23, 31, Aib32]-Pancreatic Polypeptide (Human)

Catalog #	Source	Reactivity	Applications
CCP1778	Synthetic		
Description	Peptide to [Gly1, Ser3, 22, Gln4, 34, Thr6, Ala19, Tyr21, Ala23, 31, Aib32]-Pancreatic Polypeptide (Human)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C183H281N57O54S2		
Molecular Weight	4207.73		
Purity	> 95%		
Chemical Structure	Gly - Pro - Ser - Gln - Pro - Thr - Tyr - Pro - Gly - Asp - Asn - Ala - Thr - Pro - Glu - Gln - Met - Ala - Arg - Tyr - Tyr - Ser - Ala - Leu - Arg - Arg - Tyr - Ile - Asn - Met - Ala - Aib - Arg - Gln - Arg - Tyr - NH2		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com