


**[Gly14]-Humanin (Human)**

Catalog #	Source	Reactivity	Applications
CCP1780	Synthetic		

<b>Description</b>	Peptide to [Gly14]-Humanin (Human)
<b>Form</b>	Lyophilized powder
<b>CAS Number</b>	330936-70-4
<b>Molecular Formula</b>	C118H202N34O31S2
<b>Molecular Weight</b>	2657.25
<b>Purity</b>	> 95%
<b>Chemical Structure</b>	Met - Ala - Pro - Arg - Gly - Phe - Ser - Cys - Leu - Leu - Leu - Leu - Thr - Gly - Glu - Ile - Asp - Leu - Pro - Val - Lys - Arg - Arg - Ala
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.

**Structure**



**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

**COHESION BIOSCIENCES LIMITED**

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)