

Product Data Sheet

Adrenomedullin (11-50) (Rat)

Catalog #	Source	Reactivity	Applications
CCP1870	Synthetic		
Description	Peptide to Adrenomedullin (11-50) (Rat)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C194H304N58O59S4		
Molecular Weight	4521.20		
Purity	> 95%		
Chemical Structure	Ser - Thr - Gly - Cys - Arg - Phe - Gly - Thr - Cys - Thr - Met - Gln - Lys - Leu - Ala - His - Gln - Ile - Tyr - Gln - Phe - Thr - Asp - Lys - Asp - Lys - Asp - Gly - Met - Ala - Pro - Arg - Asn - Lys - Ile - Ser - Pro - Gln - Gly - Tyr - NH2 (Disulfide)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com