

Product Data Sheet

Amyloid Bri Protein (1-34)

Catalog #	Source	Reactivity	Applications
CCP1923	Synthetic		
Description	Peptide to Amyloid Bri Protein (1-34)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C173H273N49O52S2		
Molecular Weight	3935.55		
Purity	> 95%		
Chemical Structure	Pyr - Ala - Ser - Asn - Cys - Phe - Ala - Ile - Arg - His - Phe - Glu - Asn - Lys - Phe - Ala - Val - Glu - Thr - Leu - Ile - Cys - Ser - Arg - Thr - Val - Lys - Lys - Asn - Ile - Ile - Glu - Glu - Asn (Disulfide bridge Cys5 - Cys22)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com