

Product Data Sheet

ANF (1-24) (Frog)

Catalog #	Source	Reactivity	Applications
CCP1930	Synthetic		
Description	Peptide to ANF (1-24) (Frog)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C103H165N37O34S3		
Molecular Weight	2561.87		
Purity	> 95%		
Chemical Structure	Ser - Ser - Asp - Cys - Phe - Gly - Ser - Arg - Ile - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Met - Gly - Cys - Gly - Arg - Arg - Phe (Disulfide bridge Cys4 - Cys20)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com