


Calpain Inhibitor Peptide

Catalog #	Source	Reactivity	Applications
CCP2160	Synthetic		
Description	Peptide to Calpain Inhibitor Peptide		
Form	Lyophilized powder		
CAS Number	128578-18-7		
Molecular Formula	C ₁₄₀ H ₂₂₇ N ₃₅ O ₄₄ S		
Molecular Weight	3136.64		
Purity	> 95%		
Chemical Structure	Asp - Pro - Met - Ser - Ser - Thr - Tyr - Ile - Glu - Glu - Leu - Gly - Lys - Arg - Glu - Val - Thr - Ile - Pro - Pro - Lys - Tyr - Arg - Glu - Leu - Leu - Ala		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com