

Product Data Sheet

Decorsin (Leech)

Catalog #	Source	Reactivity	Applications
CCP2244	Synthetic		
Description	Peptide to Decorsin (Leech)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C179H277N55O62S6		
Molecular Weight	4383.87		
Purity	> 95%		
Chemical Structure	Ala - Pro - Arg - Leu - Pro - Gln - Cys - Gln - Gly - Asp - Asp - Gln - Glu - Lys - Cys - Leu - Cys - Asn - Lys - Asp - Glu - Cys - Pro - Pro - Gly - Gln - Cys - Arg - Phe - Pro - Arg - Gly - Asp - Ala - Asp - Pro - Tyr - Cys - Glu		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com