


Hepatitis B Virus Receptor Binding Fragment

Catalog #	Source	Reactivity	Applications
CCP2437	Synthetic		
Description	Peptide to Hepatitis B Virus Receptor Binding Fragment		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C140H185N35O42		
Molecular Weight	3030.20		
Purity	> 95%		
Chemical Structure	Pro - Leu - Gly - Phe - Phe - Pro - Asp - His - Gln - Leu - Asp - Pro - Ala - Phe - Gly - Ala - Asn - Ser - Asn - Asn - Pro - Asp - Trp - Asp - Phe - Asn - Pro		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com