

Product Data Sheet

HIV-gp120-41-C

Catalog #	Source	Reactivity	Applications
CCP2463	Synthetic		
Description		Peptide to HIV-gp120-41-C	
Form		Lyophilized powder	
CAS Number			
Molecular Formula		C116H164N32O31S	
Molecular Weight		2534.86	
Purity		> 95%	
Chemical Structure		Val - Pro - Trp - Asn - Ser - Ser - Trp - Ser - Asn - Lys - Leu - Asp - Arg - Ile - Trp - Asn - Asn - Met - Thr - Trp	
Storage/Stability		Shipped at 4°C. Store at -20°C for one year.	
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com