


HIV-gp120-Fragment (254-274)

Catalog #	Source	Reactivity	Applications
CCP2467	Synthetic		
Description	Peptide to HIV-gp120-Fragment (254-274)		
Form	Lyophilized powder		
CAS Number	113944-46-0		
Molecular Formula	C ₉₅ H ₁₆₂ N ₂₈ O ₃₀ S		
Molecular Weight	2208.58		
Purity	> 95%		
Chemical Structure	Cys - Thr - His - Gly - Ile - Arg - Pro - Val - Val - Ser - Thr - Gln - Leu - Leu - Leu - Asn - Gly - Ser - Leu - Ala - Glu		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com