

Insulin Receptor (1142-1153) Amide

Catalog #	Source	Reactivity	Applications
CCP2525	Synthetic		
Description	Peptide to Insulin Receptor (1142-1153) Amide		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C72H108N20O23		
Molecular Weight	1621.78		
Purity	> 95%		
Chemical Structure	Thr - Arg - Asp - Ile - Tyr - Glu - Thr - Asp - Tyr - Tyr - Arg - Lys - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED**WEB**www.cohesionbio.com**ORDER**order@cohesionbio.com**SUPPORT**techsupport@cohesionbio.com**CUSTOM**custom@cohesionbio.com