

Product Data Sheet

Lytic Peptide SB-37

Catalog #	Source	Reactivity	Applications
CCP2593	Synthetic		
Description	Peptide to Lytic Peptide SB-37		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C188H320N54O45S		
Molecular Weight	4089.05		
Purity	> 95%		
Chemical Structure	Met - Pro - Lys - Trp - Lys - Val - Phe - Lys - Lys - Ile - Glu - Lys - Val - Gly - Arg - Asn - Ile - Arg - Asn - Gly - Ile - Val - Lys - Ala - Gly - Pro - Ala - Ile - Ala - Val - Leu - Gly - Glu - Ala - Lys - Ala - Leu - Gly		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com