

Product Data Sheet

Neuromedin S (Rat)

Catalog #	Source	Reactivity	Applications
CCP2660	Synthetic		
Description	Peptide to Neuromedin S (Rat)		
Form	Lyophilized powder		
CAS Number	843782-19-4		
Molecular Formula	C193H307N57O49S		
Molecular Weight	4241.97		
Purity	> 95%		
Chemical Structure	Leu - Pro - Arg - Leu - Leu - His - Thr - Asp - Ser - Arg - Met - Ala - Thr - Ile - Asp - Phe - Pro - Lys - Lys - Asp - Pro - Thr - Thr - Ser - Leu - Gly - Arg - Pro - Phe - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com