

Product Data Sheet

Neuropeptide W-30 (Human)

| Catalog # | Source | Reactivity | Applications |
|--------------------|---|------------|--------------|
| CCP2667 | Synthetic | | |
| Description | Peptide to Neuropeptide W-30 (Human) | | |
| Form | Lyophilized powder | | |
| CAS Number | 383415-80-3 | | |
| Molecular Formula | C165H249N49O37S | | |
| Molecular Weight | 3543.17 | | |
| Purity | > 95% | | |
| Chemical Structure | Trp - Tyr - Lys - His - Val - Ala - Ser - Pro - Arg - Tyr - His - Thr - Val - Gly - Arg - Ala - Ala - Gly - Leu - Leu - Met - Gly - Leu - Arg - Arg - Ser - Pro - Tyr - Leu - Trp | | |
| Storage/Stability | Shipped at 4°C. Store at -20°C for one year. | | |
| Structure | | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com