

Product Data Sheet

Neuropeptide W-30 (Rat)

Catalog #	Source	Reactivity	Applications
CCP2668	Synthetic		
Description	Peptide to Neuropeptide W-30 (Rat)		
Form	Lyophilized powder		
CAS Number	383415-90-5		
Molecular Formula	C165H249N49O38S		
Molecular Weight	3559.17		
Purity	> 95%		
Chemical Structure	Trp - Tyr - Lys - His - Val - Ala - Ser - Pro - Arg - Tyr - His - Thr - Val - Gly - Arg - Ala - Ser - Gly - Leu - Leu - Met - Gly - Leu - Arg - Arg - Ser - Pro - Tyr - Leu - Trp		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com