

Product Data Sheet

Ornipressin Acetate

Catalog #	Source	Reactivity	Applications
CCP2688	Synthetic		
Description		Peptide to Ornipressin Acetate	
Form		Lyophilized powder	
CAS Number		3397-23-7	
Molecular Formula		C45H63N13O12S2	
Molecular Weight		1042.22	
Purity		> 95%	
Chemical Structure		c[Cys - Tyr - Phe - Gln - Asn - Cys] - Pro - Orn - Gly - NH2	
Storage/Stability		Shipped at 4°C. Store at -20°C for one year.	
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com