

## Product Data Sheet

## V5 Epitope Tag

Catalog #	Source	Reactivity	Applications
CCP2974	Synthetic		
<b>Description</b>		Peptide to V5 Epitope Tag	
<b>Form</b>		Lyophilized powder	
<b>CAS Number</b>			
<b>Molecular Formula</b>		C64H1081N16O20	
<b>Molecular Weight</b>		1421.67	
<b>Purity</b>		> 95%	
<b>Chemical Structure</b>		Gly - Lys - Pro - Ile - Pro - Asn - Pro - Leu - Leu - Gly - Leu - Asp - Ser - Thr	
<b>Storage/Stability</b>		Shipped at 4°C. Store at -20°C for one year.	
<b>Structure</b>			

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)