

# Product Data Sheet

## Goat Anti-Alpaca IgG - Biotin

| Catalog #                | Source                                                                                                                                                    | Reactivity | Applications  |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| CSA2283                  | Goat                                                                                                                                                      | AI         | E, WB, IH, IC |
| <b>Description</b>       | Goat Polyclonal Secondary Antibody to Alpaca IgG Biotin labeled                                                                                           |            |               |
| <b>Immunogen</b>         | Alpaca IgG                                                                                                                                                |            |               |
| <b>Purification</b>      | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                         |            |               |
| <b>Specificity</b>       | By immunoelectrophoresis and ELISA this antibody reacts specifically with Alpaca IgG. No antibody was detected against non immunoglobulin serum proteins. |            |               |
| <b>Clonality</b>         | Polyclonal                                                                                                                                                |            |               |
| <b>Conjugation</b>       | Biotin                                                                                                                                                    |            |               |
| <b>Form</b>              | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide                                           |            |               |
| <b>Dilution</b>          | E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)                                                                     |            |               |
| <b>Storage/Stability</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                          |            |               |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)