

## Product Data Sheet

## Anti-FLAG-tag Antibody - FITC

Catalog #	Source	Reactivity	Applications
CSA5714	Rabbit	N/A	WB, IF, FC
<b>Description</b>	Rabbit polyclonal antibody to FLAG-tag FITC labeled		
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence of FLAG-tag.		
<b>Purification</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.		
<b>Specificity</b>	Recognizes C-terminal, internal, and N-terminal FLAG-tag fusion proteins.		
<b>Clonality</b>	Polyclonal		
<b>Conjugation</b>	FITC		
<b>Form</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide		
<b>Dilution</b>	WB (1/2000 - 1/5000), IF (1/500 - 1/2000), FC (1/1000 - 1/4000)		
<b>Storage/Stability</b>	Shipped at 4° C. Upon delivery aliquot and store at -20° C for one year. Avoid freeze/thaw cycles.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

WEB  
[www.cohesionbio.com](http://www.cohesionbio.com)

ORDER  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

SUPPORT  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

CUSTOM  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)