

Product Data Sheet

Streptavidin-HRP labeled

Catalog #	Source	Reactivity	Applications
CSA6215	Native	N/A	E, WB, IH, IC
Description	Streptavidin to HRP labeled		
Immunogen			
Purification			
Specificity	2-7 functional groups on each Streptavidin molecule		
Clonality			
Conjugation	HRP		
Form	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.02% Sodium Azide.		
Dilution	E (1/5000 - 1/20000), WB (1/5000 - 1/10000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com