

Anti-Beta-tubulin Antibody - AcalephFluor750

Catalog #	Source	Reactivity	Applications
CSA6512	Mouse	H, M, R	IF, FC
Description	Mouse monoclonal antibody to Beta-tubulin AcalephFluor750 labeled		
Immunogen	Recombinant protein corresponding to human Beta-tubulin.		
Purification	Affinity chromatography		
Specificity	Recognizes endogenous levels of Beta-tubulin protein.		
Clonality	Monoclonal		
Conjugation	AcalephFluor750		
Form	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide		
Dilution	IF (1/500 - 1/2000), FC (1/1000 - 1/4000)		
Gene Symbol	TUBB		
Alternative Names	TUBB5; Tubulin beta chain; Tubulin beta-5 chain		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
203068	22154	29214	
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
P07437	P99024	P69897	
Fluorophore	Ex/Em = 750/777 nm		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com