

Product Data Sheet

Sorbitol Dehydrogenase

Catalog #	Source	Reactivity	Applications
CZM1141	E.coli		
Description	Recombinant Sorbitol Dehydrogenase from E.coli		
Specificity	>40 U/mg protein, >10 U/mg powder		
Form	Light tan white powder		
Alternative Names	EC 1.1.1.14; SDH; L-iditol 2-dehydrogenase; polyol dehydrogenase; sorbitol dehydrogenase; L-iditol:NAD ⁺ 5-oxidoreductase; L-iditol (sorbitol) dehydrogenase; glucitol dehydrogenase; L-iditol:NAD ⁺ oxidoreductase; NAD ⁺ -dependent sorbitol dehydrogenase; NAD ⁺ -		
CAS Number	9028-21-1		
Molecular Weight			
Purity	>90%		
Unit Definition	One unit is defined as the amount of enzyme will reduce 1 micromole D-fructose per minute at pH 7.6 and 37°C.		
Storage/Stability	Shipped at 4°C. Store at -20°C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com