

Anti-Malaria-Pf-HRP2 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1123-C	Mouse	Malaria	SE, LFIA, CLIA
Description	Mouse monoclonal capture antibody to Malaria-Pf-HRP2 Antibody		
Immunogen	Recombinant protein of Malaria-Pf-HRP2		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant Malaria-Pf-HRP2 protein		
Clonality	Monoclonal		
Conjugation			
Form	Liquid in in PBS, pH 7.4, 0.05% Proclin 300.		
Dilution			
Gene Symbol			
Alternative Names			
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com