

Anti-RPL30 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1145-C	Mouse	H	SE
Description	Mouse monoclonal capture antibody to RPL30		
Immunogen	Recombinant protein of human RPL30		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human RPL30 protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	RPL30		
Alternative Names	60S ribosomal protein L30		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
6156			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
P62888			
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1145-D; Recommmended standard protein: CST1145		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com