

## Anti-BCL6 Antibody (Capture)

| Catalog #           | Source                                                                                                                                                                                 | Reactivity        | Applications |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|--------------|
| CAP1172-C           | Mouse                                                                                                                                                                                  | H                 | SE           |
| Description         | Mouse monoclonal capture antibody to BCL6                                                                                                                                              |                   |              |
| Immunogen           | Recombinant protein of human BCL6                                                                                                                                                      |                   |              |
| Purification        | The antibody was purified by Protein A/G affinity chromatography.                                                                                                                      |                   |              |
| Specificity         | Recognizes native and recombinant human BCL6 protein                                                                                                                                   |                   |              |
| Clonality           | Monoclonal                                                                                                                                                                             |                   |              |
| Conjugation         |                                                                                                                                                                                        |                   |              |
| Form                | Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.                                                                                                     |                   |              |
| Dilution            |                                                                                                                                                                                        |                   |              |
| Gene Symbol         | BCL6                                                                                                                                                                                   |                   |              |
| Alternative Names   | BCL5; LAZ3; ZBTB27; ZNF51; B-cell lymphoma 6 protein; BCL-6; B-cell lymphoma 5 protein; BCL-5; Protein LAZ-3; Zinc finger and BTB domain-containing protein 27; Zinc finger protein 51 |                   |              |
| Entrez Gene (Human) | Entrez Gene (Mouse)                                                                                                                                                                    | Entrez Gene (Rat) |              |
| 604                 |                                                                                                                                                                                        |                   |              |
| SwissProt (Human)   | SwissProt (Mouse)                                                                                                                                                                      | SwissProt (Rat)   |              |
| P41182              |                                                                                                                                                                                        |                   |              |
| Directions for Use  | Recommmended pair for sandwich immunoassay - Detection: CAP1172-D;<br>Recommmended standard protein: CST1172                                                                           |                   |              |
| Storage/Stability   | Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.                                                                                         |                   |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)