

Anti-SNRNP27 Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1176-C	Mouse	H	SE
Description	Mouse monoclonal capture antibody to SNRNP27		
Immunogen	Recombinant protein of human SNRNP27		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human SNRNP27 protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	SNRNP27		
Alternative Names	U4/U6.U5 small nuclear ribonucleoprotein 27 kDa protein; U4/U6.U5 snRNP 27 kDa protein; U4/U6.U5-27K; Nucleic acid-binding protein RY-1; U4/U6.U5 tri-snRNP-associated 27 kDa protein; 27K; U4/U6.U5 tri-snRNP-associated protein 3		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
11017			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q8WVK2			
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1176-D; Recommmended standard protein: CST1176		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com