

Anti-KIAA1755 Antibody (Capture)

Catalog #	Source	Reactivity	Applications		
CAP1227-C	Mouse	H	SE		
Description	Mouse monoclonal capture antibody to KIAA1755				
Immunogen	Recombinant protein of human KIAA1755				
Purification	The antibody was purified by Protein A/G affinity chromatography.				
Specificity	Recognizes native and recombinant human KIAA1755 protein				
Clonality	Monoclonal				
Conjugation					
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.				
Dilution					
Gene Symbol	KIAA1755				
Alternative Names	Uncharacterized protein KIAA1755				
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)			
85449					
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)			
Q5JYT7					
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1227-D; Recommmended standard protein: CST1227				
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.				

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com