

Anti-FAM81A Antibody (Detection)

Catalog #	Source	Reactivity	Applications		
CAP1243-D	Mouse	H	SE		
Description	Mouse monoclonal detection antibody to FAM81A				
Immunogen	Recombinant protein of human FAM81A				
Purification	The antibody was purified by Protein A/G affinity chromatography.				
Specificity	Recognizes native and recombinant human FAM81A protein				
Clonality	Monoclonal				
Conjugation					
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.				
Dilution					
Gene Symbol	FAM81A				
Alternative Names	Protein FAM81A				
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)			
145773					
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)			
Q8TBF8					
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1243-C; Recommmended standard protein: CST1243				
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.				

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com