

Anti-CIBAR1 Antibody (Capture)

Catalog #	Source	Reactivity	Applications		
CAP1335-C	Mouse	H	SE		
Description	Mouse monoclonal capture antibody to CIBAR1				
Immunogen	Recombinant protein of human CIBAR1				
Purification	The antibody was purified by Protein A/G affinity chromatography.				
Specificity	Recognizes native and recombinant human CIBAR1 protein				
Clonality	Monoclonal				
Conjugation					
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.				
Dilution					
Gene Symbol	FAM92A1				
Alternative Names	Protein FAM92A1				
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)			
137392					
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)			
A1XBS5					
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1335-D; Recommmended standard protein: CST1335				
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.				

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com