

Anti-USP6NL Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1337-C	Mouse	H	SE
Description	Mouse monoclonal capture antibody to USP6NL		
Immunogen	Recombinant protein of human USP6NL		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human USP6NL protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	USP6NL		
Alternative Names	KIAA0019; USP6 N-terminal-like protein; Related to the N-terminus of tre; RN-tre		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
9712			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q92738			
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1337-D; Recommmended standard protein: CST1337		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com