

Anti-ZNF788P Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1354-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to ZNF788P		
Immunogen	Recombinant protein of human ZNF788P		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human ZNF788P protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	ZNF788		
Alternative Names	Zinc finger protein 788		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
388507			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q6ZQV5			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1354-C; Recommmended standard protein: CST1354		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com