

Anti-HYCC2 Antibody (Detection)

Catalog #	Source	Reactivity	Applications		
CAP1364-D	Mouse	H	SE		
Description	Mouse monoclonal detection antibody to HYCC2				
Immunogen	Recombinant protein of human HYCC2				
Purification	The antibody was purified by Protein A/G affinity chromatography.				
Specificity	Recognizes native and recombinant human HYCC2 protein				
Clonality	Monoclonal				
Conjugation					
Form	Mouse IgG3. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.				
Dilution					
Gene Symbol	FAM126B				
Alternative Names	Protein FAM126B				
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)			
285172					
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)			
Q8IXS8					
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1364-C; Recommmended standard protein: CST1364				
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.				

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com