

Anti-INTS6L Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1387-C	Mouse	H	SE
Description	Mouse monoclonal capture antibody to INTS6L		
Immunogen	Recombinant protein of human INTS6L		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human INTS6L protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	DDX26B		
Alternative Names	Protein DDX26B		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
203522			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q5JSJ4			
Directions for Use	Recommmended pair for sandwich immunoassay - Detection: CAP1387-D; Recommmended standard protein: CST1387		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com