

## Anti-USP54 Antibody (Detection)

| Catalog #           | Source                                                                                                 | Reactivity        | Applications |
|---------------------|--------------------------------------------------------------------------------------------------------|-------------------|--------------|
| CAP1389-D           | Mouse                                                                                                  | H                 | SE           |
| Description         | Mouse monoclonal detection antibody to USP54                                                           |                   |              |
| Immunogen           | Recombinant protein of human USP54                                                                     |                   |              |
| Purification        | The antibody was purified by Protein A/G affinity chromatography.                                      |                   |              |
| Specificity         | Recognizes native and recombinant human USP54 protein                                                  |                   |              |
| Clonality           | Monoclonal                                                                                             |                   |              |
| Conjugation         |                                                                                                        |                   |              |
| Form                | Mouse IgG2b. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.                    |                   |              |
| Dilution            |                                                                                                        |                   |              |
| Gene Symbol         | USP54                                                                                                  |                   |              |
| Alternative Names   | C10orf29; Inactive ubiquitin carboxyl-terminal hydrolase 54; Inactive ubiquitin-specific peptidase 54  |                   |              |
| Entrez Gene (Human) | Entrez Gene (Mouse)                                                                                    | Entrez Gene (Rat) |              |
| 159195              |                                                                                                        |                   |              |
| SwissProt (Human)   | SwissProt (Mouse)                                                                                      | SwissProt (Rat)   |              |
| Q70EL1              |                                                                                                        |                   |              |
| Directions for Use  | Recommmended pair for sandwich immunoassay - Capture: CAP1389-C; Recommmeded standard protein: CST1389 |                   |              |
| Storage/Stability   | Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.         |                   |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)