

Anti-PLCH2 Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1407-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to PLCH2		
Immunogen	Recombinant protein of human PLCH2		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human PLCH2 protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2b. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	PLCH2		
Alternative Names	KIAA0450; PLCL4; 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase eta-2; Phosphoinositide phospholipase C-eta-2; Phosphoinositide phospholipase C-like 4; PLC-L4; Phospholipase C-like protein 4; Phospholipase C-eta-2; PLC-eta2		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
9651			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
O75038			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1407-C; Recommmeded standard protein: CST1407		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com