

## Anti-MTFR1 Antibody (Capture)

| Catalog #           | Source                                                                                                       | Reactivity        | Applications |
|---------------------|--------------------------------------------------------------------------------------------------------------|-------------------|--------------|
| CAP1462-C           | Mouse                                                                                                        | H                 | SE           |
| Description         | Mouse monoclonal capture antibody to MTFR1                                                                   |                   |              |
| Immunogen           | Recombinant protein of human MTFR1                                                                           |                   |              |
| Purification        | The antibody was purified by Protein A/G affinity chromatography.                                            |                   |              |
| Specificity         | Recognizes native and recombinant human MTFR1 protein                                                        |                   |              |
| Clonality           | Monoclonal                                                                                                   |                   |              |
| Conjugation         |                                                                                                              |                   |              |
| Form                | Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.                           |                   |              |
| Dilution            |                                                                                                              |                   |              |
| Gene Symbol         | MTFR1                                                                                                        |                   |              |
| Alternative Names   | CHPPR; FAM54A2; KIAA0009; Mitochondrial fission regulator 1; Chondrocyte protein with a poly-proline region  |                   |              |
| Entrez Gene (Human) | Entrez Gene (Mouse)                                                                                          | Entrez Gene (Rat) |              |
| 9650                |                                                                                                              |                   |              |
| SwissProt (Human)   | SwissProt (Mouse)                                                                                            | SwissProt (Rat)   |              |
| Q15390              |                                                                                                              |                   |              |
| Directions for Use  | Recommmended pair for sandwich immunoassay - Detection: CAP1462-D;<br>Recommmended standard protein: CST1462 |                   |              |
| Storage/Stability   | Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.               |                   |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)