

Anti-Fibrinogen gamma Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1470-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to Fibrinogen gamma		
Immunogen	Recombinant protein of human Fibrinogen gamma		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human Fibrinogen gamma protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	FGG		
Alternative Names	Fibrinogen gamma chain		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
2266			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
P02679			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1470-C; Recommmended standard protein: CST1470		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com