

Anti-CZIB Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1501-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to CZIB		
Immunogen	Recombinant protein of human CZIB		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human CZIB protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	C1orf123		
Alternative Names	UPF0587 protein C1orf123		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
54987			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q9NWW4			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1501-C; Recommmended standard protein: CST1501		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com