

Anti-Podocin Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1606-D	Mouse	H	SE
Description	Mouse monoclonal detection antibody to Podocin		
Immunogen	Recombinant protein of human Podocin		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human Podocin protein		
Clonality	Monoclonal		
Conjugation			
Form	Mouse IgG2b. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.		
Dilution			
Gene Symbol	NPHS2		
Alternative Names	Podocin		
Entrez Gene (Human)	Entrez Gene (Mouse)	Entrez Gene (Rat)	
7827			
SwissProt (Human)	SwissProt (Mouse)	SwissProt (Rat)	
Q9NP85			
Directions for Use	Recommmended pair for sandwich immunoassay - Capture: CAP1606-C; Recommmended standard protein: CST1606		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com