

# Product Data Sheet

## Anti-CD319 Antibody (Detection)

| Catalog #           | Source                                                                                                                                                | Reactivity        | Applications |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|--------------|
| CAP2066-D           | Mouse                                                                                                                                                 | H                 | SE           |
| Description         | Mouse monoclonal detection antibody to CD319                                                                                                          |                   |              |
| Immunogen           | Recombinant protein of human CD319                                                                                                                    |                   |              |
| Purification        | The antibody was purified by Protein A/G affinity chromatography.                                                                                     |                   |              |
| Specificity         | Recognizes native and recombinant human CD319 protein                                                                                                 |                   |              |
| Clonality           | Monoclonal                                                                                                                                            |                   |              |
| Conjugation         |                                                                                                                                                       |                   |              |
| Form                | Mouse IgG2a. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide.                                                                   |                   |              |
| Dilution            |                                                                                                                                                       |                   |              |
| Gene Symbol         | SLAMF7                                                                                                                                                |                   |              |
| Alternative Names   | CS1; SLAM family member 7; CD2 subset 1; CD2-like receptor-activating cytotoxic cells; CRACC; Membrane protein FOAP-12; Novel Ly9; Protein 19A; CD319 |                   |              |
| Entrez Gene (Human) | Entrez Gene (Mouse)                                                                                                                                   | Entrez Gene (Rat) |              |
| 57823               |                                                                                                                                                       |                   |              |
| SwissProt (Human)   | SwissProt (Mouse)                                                                                                                                     | SwissProt (Rat)   |              |
| Q9NQ25              |                                                                                                                                                       |                   |              |
| Directions for Use  | Recommmended pair for sandwich immunoassay - Capture: CAP2066-C;<br>Recommmended standard protein: CST2066                                            |                   |              |
| Storage/Stability   | Shipped at 4°C. Upon delivery aliquot and store at -20°C for 1 year. Avoid freeze/thaw cycles.                                                        |                   |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

### WEB

[www.cohesionbio.com](http://www.cohesionbio.com)

### ORDER

[order@cohesionbio.com](mailto:order@cohesionbio.com)

### SUPPORT

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

### CUSTOM

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)